



## IHE Integration Statement



## Aixplorer® MACH30 and MACH20 Ultrasound System V2.X software

**Revision:** A  
**Date:** August 19, 2019

**Doc Ref:** RD.DD.359 (RevA)  
**Software Version:** V2.x

## Introduction

Integrating the Healthcare Enterprise® (IHE) is an initiative by healthcare professionals and industry to improve the way computer systems in healthcare share information. Systems that support IHE Integration Profiles work together better, are easier to implement, and help care providers use information more effectively. The goal is efficient delivery of optimal patient care.

IHE Integration Statements are documents prepared and published by vendors to describe the intended conformance of their products with the IHE Technical Framework. Integration Statements tell consumers the specific IHE capabilities a given product is designed to support in terms of the key concepts of IHE: Actors and Integration Profiles.

## **IHE Integration Statement**

<b>IHE Integration Statement</b>			
Vendor	Product Name	Date	August 19, 2019
SuperSonic Imagine	Aixplorer Mach30	Version	V2.x
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles	Actors	Options	
Radiology domain			
Scheduled Workflow Scheduled Workflow.B	Acquisition Modality	Broad Worklist Query Patient Based Worklist Query PPS Exception Management	
Patient Information Reconciliation	Acquisition Modality		
Evidence Documents	Acquisition Modality		
Portable Data for Imaging	Portable Media Creator	Web Content	
Cardiology domain			
Echocardiography Workflow	Acquisition Modality	Broad Worklist Query Patient Based Worklist Query PPS Exception Management	
Evidence Documents	Acquisition Modality	Echo Evidence	
Image-Enabled Office	Acquisition Modality	Broad Worklist Query Patient Based Worklist Query PPS Exception Management	
IT infrastructure domain			
Consistent Time	Time Client		
Internet address for vendor's IHE information: <a href="http://www.supersonicimagine.com/ihe">http://www.supersonicimagine.com/ihe</a>			
Links to Standards Conformance Statements for the Implementation			
HL7	Not Applicable		
DICOM	<a href="http://www.supersonicimagine.com/dicom">http://www.supersonicimagine.com/dicom</a>		
Links to general information on IHE			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">www.ihe-j.org</a>	