



The Theragnostic Company™

## IHE Integration Statement



## Aixplorer® Ultrasound System

**Revision:** 01

**Date:** January 17 2012

**Doc Ref:** SSID02514-01

**Software Version:** V4.3

<b>IHE Integration Statement</b>			
Vendor	Product Name	Version	Date
SuperSonic Imagine	Aixplorer	V4.3	Jan 2012
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Integration Profile Option Implemented	
Scheduled Workflow	Acquisition Modality	Broad Worklist Query PPS Exception Management	
Patient Information Reconciliation	Acquisition Modality		
Internet address for vendors IHE information: <a href="http://www.supersonicimagine.fr/ihe">http://www.supersonicimagine.fr/ihe</a>			
Links to Standards Conformance Statements for the Implementation			
HL7			
DICOM	<a href="http://www.supersonicimagine.fr/dicom">http://www.supersonicimagine.fr/dicom</a>		
Links to general information on IHE			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	